

Input DNA



Desire:

1. Clean
2. Good yield
3. High molecular weight (HMW)



Why?

Assessments?

Input DNA

Bacteria:

Precipitation-based kits (Promega Wizard kit)
HMW extraction kits (Monarch, others)
column methods shear DNA, not good



Eukarya:

gentle/old school methods
Phenol/Chloroform
Shepherd's crook
Qiagen genomic tip, HMW kits




post extraction size selection:

Needle shear ultra-long DNA first
then 0.4X SPRI bead clean up

or Blue Pippin size selection
(generally not needed for Bacteria)

Clean DNA



Why clean DNA?

1. Chemicals impact downstream manipulations (library prep)
2. Chemicals or co-precipitated protein may impact proteins on flowcell

Nanopores are proteins

Need to balance amount and size

Why?

Goal is to have pores constantly sequencing DNA


Small DNA = preferentially sequenced, but better overall yield

Too (HMW) DNA = ultralong DNA clogs pores quickly

underloaded/little DNA = pores open longer, degrade faster

If HMW DNA, can fragment

Too much DNA = generally not a big problem, can help with yield



No "one size fits all solution"

Depends on genome size and amount of repeats

Depends on use of genome sequence

But

Can always dilute sample or make HMW smaller, but not the opposite


Methods for assessing DNA

Nanodrop




$260/280 \geq 1.8$
 $260/230 \geq 2.0$

Gel



~Size, ~yield, ~quality

Qubit



Yield

Bioanalyzer



Size

Flongle



~Size, ~yield, ~quality

Library preparation: ligation kit



- + : generally highest yield (>100-150 Gb on PromethION)
- : need other kits (end-repair & nick repair); cost; time
- +/- : generally high seq length, but does not get ultralong reads

DATA OUTPUT

Estimated bases

129.06 Gb

Reads generated

12.59 M

Estimated N50

20.59 kb

Total data produced (pass / fail)

1.6 TB

BASECALLING

Reads called

96.7%

Bases called (min Q score: 10)

97.06 Gb

24.46 Gb

Pass

Fail

Legend


■ Basecalled — Estimated

Estimated N50

20.59 kb

% Basecalled

96.7%



Library preparation: rapid kit



- + : faster, easier, and cheaper
- : generally lower yield (~50-80 GB on P2);
- : less control in fragmenting DNA, not good for amplicons if whole length is needed
- +/- : generally long reads, can get ultralong reads

DATA OUTPUT

Estimated bases

34.57 Gb

Reads generated

7.8 M

Estimated N50

11.8 kb

Total data produced (pass / fail)

423.54 GB

BASECALLING

Reads called

100%

Bases called (min Q score: 10)

17.68 Gb

Pass

5.18 Gb

Fail

Legend


■ Basecalled — Estimated

Estimated N50


11.8 kb

% Basecalled

100%




Overview of ligation kit




~9 hrs

Overview of rapid kit




B Simplified transposon



C End sequences

OE: CTG**A**CTCTT**A**TACACA**A**GT
 IE: CTG**T**CTCTT**G**ATCAG**A**TCT
 ME: CTG**T**CTCTT**A**TACACA**T**CT



~2 hrs

Quality control

What is the end goal?

1. Infer taxonomic identity
2. Functional genomics
3. Genome evolution (structure & variation)
4. MGEs
5. Population genomics

What do you know about the genome(s) of interest?


1. Size
2. Complexity (repeat number and size)
3. Segmental duplications
4. Ploidy

Can't do much w/ bad reads


Findings based on genome analyses should be framed on the quality of assemblies

Three important properties of reads

1. Read length




2. Accuracy



3. Evenness


Average



Low



High



Tools to generate summary statistics

QUAST: quality assessment tool for genome assemblies

Alexey Gurevich^{1,*}, Vladislav Saveliev¹, Nikolay Vyahhi¹ and Glenn Tesler²

MultiQC: summarize analysis results for multiple tools and samples in a single report

Philip Ewels^{1,*}, Måns Magnusson², Sverker Lundin³ and Max Käller³

QUAST: helps assess quality of genome assembly

MultiQC: gathers, summarizes, and visualizes stats on reads &/or assemblies

Equipment and computing needs

Magnetic rack for bead cleanup (cheapest to 3D print or order from (Sergi Lab Supply)

<https://3d.nih.gov/entries/3DPX-013579>

<https://www.kjmagnetics.com/proddetail.asp?prod=BZZ082>

Qubit not required but highly recommended

Linux computer with GPU for running the sequencer

GPU is the most important part, needed for basecalling

NVIDIA RTX 4090 (or better) is good enough for P2 solo

System76 makes good computers with Linux pre-installed

ONT all-in-one devices (MinIT, Mk1c, P2i...) are convenient but tend to get outdated quickly...

Magnetic rack for microtubes

Za ZaneColaric

